


Product Datasheet

MGST3 antibody (orb584646)

Description	Rabbit polyclonal antibody to MGST3
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	The immunogen is a synthetic peptide directed towards the N terminal region of human MGST3
Target	MGST3
Preservatives	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration	0.5 mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Protein Sequence	Synthetic peptide located within the following region: LAINVSKARKKYYKVEYPIMYSTDPENGHIFNCIQRAHQNTLEVYPPFLFF
Clonality	Polyclonal
MW	16kDa
Uniprot ID	O14880
NCBI	NP_004519
Expiration Date	12 months from date of receipt.



Rabbit Anti-MGST3 antibody, Formalin Fix...



WB Suggested Anti-MGST3 Antibody Titration...