

Product Datasheet

CENP-A (phospho-Ser7) antibody (orb5847)

Description	CENP-A (phospho-Ser7) antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC-Fr, WB
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human CENP-A around the phosphorylation site of Ser7 (RR(p-S)RK)
Target	Phospho-CENPA (Ser7)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	15 kDa
Uniprot ID	P49450
Expiration Date	12 months from date of receipt.