

Product Datasheet


Hspd1 antibody (orb584895)

Description	Rabbit polyclonal antibody to Hspd1
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Target	Hspd1
Preservatives	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration	0.5 mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Protein Sequence	Synthetic peptide located within the following region: LDITTSEYEKEKLNERLAKLSDGVAVLKVGGTSDVEVNEKKDRVTDALNA
Clonality	Polyclonal
MW	63kDa
Uniprot ID	P63038
NCBI	NP_071565
Expiration Date	12 months from date of receipt.




Host: Mouse,
Target Name:
HSPD1, Sample

...



Host: Rabbit,
Target Name:
HSPD1,
Sample...



Rabbit Anti-
HSPD1
Antibody,
Catalog Numb...