

Product Datasheet

INPP5K antibody (orb584915)

Description	Rabbit polyclonal antibody to INPP5K
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	The immunogen is a synthetic peptide directed towards the C terminal region of human INPP5K
Target	INPP5K
Preservatives	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration	0.5 mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Protein Sequence	Synthetic peptide located within the following region: PPTYKFDRNSNDYDTSEKRRKPAWTRILWRLKRQPCAGPDTPIPPASHF
Clonality	Polyclonal
MW	49kDa
Uniprot ID	Q9BT40
NCBI	NP_057616
Expiration Date	12 months from date of receipt.



WB Suggested
Anti-INPP5K
Antibody,
Titra...