



Product Datasheet

EIF6 antibody (orb585008)


Description	Rabbit polyclonal antibody to EIF6
Species/Host	Rabbit
Reactivity	Frog, Human
Conjugation	Unconjugated
Tested Applications	WB
Target	EIF6
Preservatives	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration	0.5 mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Protein Sequence	Synthetic peptide located within the following region: AVRASFENNCEIGCFKLTNTYCLVAIGGSENFYSVFEGELSDTIPVVHA
Clonality	Polyclonal
MW	27kDa
Uniprot ID	P56537
NCBI	NP_002203
Expiration Date	12 months from date of receipt.



Host: Rabbit,
Target Name:
EIF6, Sample ...



Rabbit Anti-EIF6
antibody,
Formalin Fixe...



Rabbit Anti-EIF6
antibody,
Formalin Fixe...