



Product Datasheet

CUL4A antibody (orb585081)


Description	Rabbit polyclonal antibody to CUL4A
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, WB
Target	CUL4A
Preservatives	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration	0.5 mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Protein Sequence	Synthetic peptide located within the following region: SLACGKARVLIKSPKGKEVEDGDKFIFNGEFKHKLFRIKINQIQMKETVE
Clonality	Polyclonal
MW	88kDa
Uniprot ID	Q13619
NCBI	NP_001008895
Expiration Date	12 months from date of receipt.



Host: Rabbit,
Target Name:
CUL4A,
Sample...



Host: Rabbit,
Target Name:
CUL4A,
Sample...



Rabbit Anti-
CUL4A
Antibody,
Catalog Numb...