



Product Datasheet

BRCC3 antibody (orb585188)


Description	Rabbit polyclonal antibody to BRCC3
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF, WB
Target	BRCC3
Preservatives	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration	0.5 mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Protein Sequence	Synthetic peptide located within the following region: HNGSVFTKNLCSQMSAVSGPLLQWLEDRLEQNQQHLQELQQEKEELMQEL
Clonality	Polyclonal
MW	35kDa
Uniprot ID	P46736
NCBI	NP_077308
Expiration Date	12 months from date of receipt.



Host: Rabbit,
Target Name:
BRCC3,
Sample...



Sample Type:
MCF7, Primary
Antibody dilu...



WB Suggested
Anti-BRCC3
Antibody,
Titrat...