



Product Datasheet

UBE2Z antibody (orb585478)


Description	Rabbit polyclonal antibody to UBE2Z
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Target	UBE2Z
Preservatives	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration	0.5 mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Protein Sequence	Synthetic peptide located within the following region: GEKRGHFDYQSLLMRLGLIRQKVLRLHNENAEMDSDSSSSGTETDLHGS
Clonality	Polyclonal
MW	27kDa
Uniprot ID	Q9H832
NCBI	NP_075567
Expiration Date	12 months from date of receipt.



Host: Rabbit,
Target Name:
UBE2Z,
Sample...



Host: Rabbit,
Target Name:
UBE2Z,
Sample...



WB Suggested
Anti-UBE2Z
Antibody,
Titrat...