



Product Datasheet

iNOS antibody (orb585494)


Description	Rabbit polyclonal antibody to iNOS
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, WB
Target	NOS2
Preservatives	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration	0.5 mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Protein Sequence	Synthetic peptide located within the following region: PCATSSPVTQDDLQYHNLSKQNESPQLVETGKKSPESLVKLDATPLSS
Clonality	Polyclonal
MW	110kDa
Uniprot ID	P35228
NCBI	EAW51048
Expiration Date	12 months from date of receipt.



Lanes: Lane 1:
30 ug human
placental tis...



Rabbit Anti-
NOS2 Antibody,
Catalog
Numbe...



Application:
Western
blotting,
Species+t...