



Product Datasheet

CD155 antibody (orb585504)

Description	Rabbit polyclonal antibody to CD155
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, WB
Target	PVR
Preservatives	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration	0.5 mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Protein Sequence	Synthetic peptide located within the following region: FHQTQGPSYSESKRLEFVAARLGAELRNASLRMFGRLRVEDEGNYTCLFVT
Clonality	Polyclonal
MW	45kDa
Uniprot ID	P15151
NCBI	NP_006496
Expiration Date	12 months from date of receipt.



Host: Rabbit,
Target: PVR,
Positive cont...



Application: IHC,
Species+tissue/cell
ty...



WB Suggested
Anti-PVR Antibody,
Titratio...