



Product Datasheet

ADIPOR1 antibody (orb585607)


Description	Rabbit polyclonal antibody to ADIPOR1
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IHC, WB
Target	ADIPOR1
Preservatives	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration	0.5 mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Protein Sequence	Synthetic peptide located within the following region: FPGKFDIWFQSHQIFHVLVAAAFVHFYGVSNLQEFYRGLGEGCTDDTLL
Clonality	Polyclonal
MW	41 kDa
Uniprot ID	Q96A54
NCBI	NP_057083
Expiration Date	12 months from date of receipt.



25 ug of the indicated Human whole cell ...



Host: Mouse,
Target Name:
ADIPOR1,
Sampl...



Rabbit Anti-ADIPOR1 Antibody, Catalog Nu...