


Product Datasheet

COQ9 antibody (orb585618)

Description	Rabbit polyclonal antibody to COQ9
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, WB
Target	COQ9
Preservatives	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration	0.5 mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Protein Sequence	Synthetic peptide located within the following region: ASAVGLRSSDEQKQPPNSFSQHQHSETQGAEKPDPESSHSPPRYTDQGGE
Clonality	Polyclonal
MW	35kDa
Uniprot ID	O75208
NCBI	NP_064708
Expiration Date	12 months from date of receipt.



Rabbit Anti-COQ9 Antibody, Catalog Numbe...



WB Suggested Anti-COQ9 Antibody, Titrati...