



Product Datasheet

IL26 antibody (orb586073)


Description	Rabbit polyclonal antibody to IL26
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Target	IL26
Preservatives	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration	0.5 mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Protein Sequence	Synthetic peptide located within the following region: LVTLSLAIAKHKQSSFTKSCYPRGTLTQAVDLYIKAAWLKATIPEDRIK
Clonality	Polyclonal
MW	20 kDa
Uniprot ID	Q9NPH9
NCBI	NP_060872
Expiration Date	12 months from date of receipt.



25 ug of the indicated Human whole cell ...



Host: Rabbit, Target Name: IL26, Sample ...



Host: Rabbit, Target Name: IL26, Sample ...