



Product Datasheet

GPM6A antibody (orb586134)


Description	Rabbit polyclonal antibody to GPM6A
Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Target	GPM6A
Preservatives	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration	0.5 mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Protein Sequence	Synthetic peptide located within the following region: IAMVHYLMVLSANWAYVKDACRMQKYEDIKSKEEQELHDIHSTRSKERLN
Clonality	Polyclonal
MW	31kDa
Uniprot ID	P51674
NCBI	NP_005268
Expiration Date	12 months from date of receipt.



Host: Rabbit,
Target Name:
GPM6A,
Antibo...



Host: Rabbit,
Target Name:
GPM6A,
Sample...



Host: Rabbit,
Target: GPM6A,
Positive co...