


Product Datasheet

ANKFY1 antibody (orb586584)

Description	Rabbit polyclonal antibody to ANKFY1
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	The immunogen for Anti-ANKFY1 antibody is: synthetic peptide directed towards the N-terminal region of Human ANFY1
Target	ANKFY1
Preservatives	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration	0.5 mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Protein Sequence	Synthetic peptide located within the following region: TPLHLVALYSSKKHSADVMSEMAQIAEALLQAGANPNMQDSKGRTPHVS
Clonality	Polyclonal
MW	128 kDa
Uniprot ID	Q9P2R3
NCBI	XP_005256736.1
Expiration Date	12 months from date of receipt.



Host: Rabbit,
Target Name:
ANKFY1,
Sample...