

Product Datasheet


FAM108A1 antibody (orb586821)

Description	Rabbit polyclonal antibody to FAM108A1
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	The immunogen is a synthetic peptide directed towards the C-terminal region of Human FAM108A1
Target	ABHD17A
Preservatives	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration	0.5 mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Protein Sequence	Synthetic peptide located within the following region: VPTVDLASRYECAAVVLHSPLTSGMRVAFPDTKKTYCFDAFPNIEKVSKI
Clonality	Polyclonal
MW	39kDa
Uniprot ID	H2R6I8
NCBI	NP_112490
Expiration Date	12 months from date of receipt.




Host: Rabbit,
Target Name:
FAM108A1,
Sam...

FAM108A1



Host: Rabbit,
Target Name:
FAM108A1,
Sam...

FAM108A1



Host: Rabbit,
Target Name:
FAM108A1,
Sam...