

Product Datasheet

PDZD11 antibody (orb586877)

Description	Rabbit polyclonal antibody to PDZD11
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	The immunogen is a synthetic peptide directed towards the C-terminal region of Human PDZD11
Target	PDZD11
Preservatives	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration	0.5 mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Protein Sequence	Synthetic peptide located within the following region: EGDQVLAVNDVDFQDIEHSCAVEILKTAREISMRVRRFFPYNYHRQKERTV
Clonality	Polyclonal
MW	15kDa
Uniprot ID	Q5EBL8
NCBI	NP_057568
Expiration Date	12 months from date of receipt.



Sample Type:
OVCAR-3 Whole
cell lysates,...