



Product Datasheet

Cytokeratin 17 (phospho-Ser44) antibody (orb5870)

Description	Cytokeratin 17 (phospho-Ser44) antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ICC, WB
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human CK 17 around the phosphorylation site of Ser44 (LG(P-S)A)
Target	Phospho-Cytokeratin 17(Ser44)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	47 kDa
Uniprot ID	Q04695
Dilution Range	WB=1:500-2000, ICC=1:100, IHC-F=1:100-500, IF=1:100-500, ELISA=1:5000-10000
Expiration Date	12 months from date of receipt.



ICC analysis of HeLa cell using Phospho-...



Western blot: Lane 1: Human HeLa cell; Lane...