



Product Datasheet

PTAR1 antibody (orb587726)

Description	Rabbit polyclonal antibody to PTAR1
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	The immunogen is a synthetic peptide directed towards the C-terminal region of Human PTAR1
Target	PTAR1
Preservatives	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration	0.5 mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Protein Sequence	Synthetic peptide located within the following region: LISQTVIDSSVMEQNPLRSEPALVPPKDEEAAVSTEEPRINLPHLLLEEV
Clonality	Polyclonal
MW	46kDa
Uniprot ID	Q7Z6K3
NCBI	NP_001093136
Expiration Date	12 months from date of receipt.



Host: Rabbit,
Target Name:
PTAR1,
Sample...



Host: Rabbit,
Target Name:
PTAR1,
Sample...