

Product Datasheet

Claudin 6 (phospho-Tyr219) antibody (orb5880)

Description	Claudin 6 (phospho-Tyr219) antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC-Fr, IHC-P, WB
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human Claudin 6 around the phosphorylation site of Tyr219 (KN(p-Y)V)
Target	phospho-Claudin 6 (Tyr219)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	23 kDa
Uniprot ID	P56747
Expiration Date	12 months from date of receipt.