


Product Datasheet

RAD51B antibody (orb589203)

Description	Rabbit polyclonal antibody to RAD51B
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	The immunogen is a synthetic peptide directed towards the middle region of human RAD51B
Target	RAD51B
Preservatives	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration	0.5 mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Protein Sequence	Synthetic peptide located within the following region: EFDAQLQGNLKERNKFLAREASSLKYLAEFFSIPVILTNTQITTHLSGALA
Clonality	Polyclonal
MW	30 kDa
Uniprot ID	A0A0A0
NCBI	NP_002868.1
Expiration Date	12 months from date of receipt.



Host: Rabbit,
Target Name:
RAD51B,
Samp...