



Product Datasheet

CHRFAM7A antibody (orb592778)

Description	Rabbit polyclonal antibody to CHRFAM7A
Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	The immunogen is a synthetic peptide directed towards the N terminal region of human CHRNA7
Target	CHRNA7
Preservatives	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration	1.0 mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Protein Sequence	Synthetic peptide located within the following region: QGEFQRKLYKELVKNYNPLRPVANDSQPLTVYFSLSLQIMDVDEKNQV
Clonality	Polyclonal
MW	56 kDa
Uniprot ID	P36544
NCBI	NP_000737
Expiration Date	12 months from date of receipt.




25 ug of the indicated Human whole cell ...



Host: Rabbit,
Target Name:
CHRNA7,
Sampl...

CHRNA7



Host: Rabbit,
Target Name:
CHRNA7,
Sampl...