


Product Datasheet

NFATC1 antibody (orb592967)


Description	Rabbit polyclonal antibody to NFATC1
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	The immunogen is a synthetic peptide directed towards the N terminal region of human NFATC1
Target	NFATC1
Preservatives	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration	0.5 mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Protein Sequence	Synthetic peptide located within the following region: MPSTSPVPVPSKFPLGPAAAVFGRGETLGPAPRAGGTMKSAEEEHYGYASS
Clonality	Polyclonal
MW	78kDa
Uniprot ID	O95644
NCBI	NP_765978
Expiration Date	12 months from date of receipt.



Rabbit Anti-NFATC1 Antibody, Catalog Num...



WB Suggested Anti-NFATC1 Antibody Titrat...



WB Suggested Anti-NFATC1 antibody Titrat...