



Product Datasheet

GRK1 antibody (orb593282)


Description	Rabbit polyclonal antibody to GRK1
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC, WB
Immunogen	Recombinant fusion protein of human GRK1(NP_002920.1).
Preservatives	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	65kDa
Uniprot ID	Q15835
Dilution Range	WB: 1:200-2000, IHC-P: 1:50-200, IF: 1:50-200
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of mouse retina...



Immunofluorescence analysis of rat retina...



Western blot analysis of extracts of Y79...