



Product Datasheet

NFkB p100 (Phospho-S866) antibody (orb593301)


Description	Rabbit polyclonal antibody to NFkB p100 (Phospho-S866)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Immunogen	A phospho specific peptide corresponding to residues surrounding S866 of human NF-κB2
Preservatives	Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purification
MW	120kDa
Uniprot ID	Q00653
Entrez	4791
Dilution Range	WB: 1:500-2000, IHC-P: 1:50-100, IF: 1:50-100
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of methanol-
...



Immunohistochemical staining of human
br...



Western blot analysis of extracts of MDA...