

Product Datasheet

PPP1R12A (Phospho-T853) antibody (orb593316)

Description	Rabbit polyclonal antibody to PPP1R12A (Phospho-T853)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A synthetic phosphorylated peptide around T853 of human PPP1R12A.
Preservatives	Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purification
MW	125kDa
Uniprot ID	O14974
Entrez	4659
Dilution Range	WB: 1:500-2000
Expiration Date	12 months from date of receipt.