

Product Datasheet

CEBPB (Phospho-T235) antibody (orb593352)

Description	Rabbit polyclonal antibody to CEBPB (Phospho-T235)
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A synthetic phosphorylated peptide around T235 of human C/EBPB(NP_005185.2).
Preservatives	PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	38KDa
Uniprot ID	P17676
Dilution Range	WB: 1:500-2000
Expiration Date	12 months from date of receipt.