


Product Datasheet

P70S6K1 (Phospho-T389) antibody (orb593356)

Description	Rabbit polyclonal antibody to P70S6K1 (Phospho-T389)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	A phospho synthetic peptide of human P70S6K1.
Preservatives	Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purification
MW	70kDa
Uniprot ID	P23443
Entrez	6198
Dilution Range	WB: 1:500-2000
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of HeL...