

Product Datasheet

Histone H2A (Phospho-S129) antibody (orb593360)

Description	Rabbit polyclonal antibody to Histone H2A (Phospho-S129)
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	IHC
Immunogen	A synthetic phosphorylated peptide around S129 of saccharomyces cerevisiae Histone H2A(NP_010511.3?).
Preservatives	PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	14kDa
Uniprot ID	P04911
Dilution Range	IHC-P: 1:50-200
Expiration Date	12 months from date of receipt.