

## Product Datasheet

# Sabin Poliovirus Type I antibody (HRP) (orb593609)

<b>Description</b>	Bovine polyclonal antibody against Sabin Poliovirus Type I conjugated to HRP
<b>Reactivity</b>	Bovine
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	ELISA
<b>Immunogen</b>	Inactivated Sabin Poliovirus Type I
<b>Preservatives</b>	0.1M PBS, PH 8.0, with protein stabilizer
<b>Storage</b>	Store at 4°C for 6 months. It is recommended to aliquote to ~10ul per vial and store at -20°C for long term storage. Avoid repeated freezing and thawing cycles.
<b>Note</b>	For research use only
<b>Expiration Date</b>	12 months from date of receipt.