

Product Datasheet

Sabin Poliovirus Type II antibody (HRP) (orb593613)

Description	Guinea pig polyclonal antibody against Sabin Poliovirus Type II conjugated to HRP
Reactivity	Guinea pig
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Inactivated Sabin Poliovirus Type II
Preservatives	0.1M PBS, PH 8.0, with protein stabilizer
Storage	Store at 4°C for 6 months. It is recommended to aliquote to ~10ul per vial and store at -20°C for long term storage. Avoid repeated freezing and thawing cycles.
Note	For research use only
Expiration Date	12 months from date of receipt.