

Product Datasheet

Sabin Poliovirus Type I antibody (HRP) (orb593615)

Description	Rabbit polyclonal antibody against Sabin Poliovirus Type I conjugated to HRP
Reactivity	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Inactivated Sabin Poliovirus Type I
Preservatives	0.1M PBS, PH 8.0, with protein stabilizer
Storage	Store at 4°C for 6 months. It is recommended to aliquote to ~10ul per vial and store at -20°C for long term storage. Avoid repeated freezing and thawing cycles.
Note	For research use only
Expiration Date	12 months from date of receipt.