

Product Datasheet

Haemophilus Influenza Type b (Hib) antibody (HRP) (orb593626)

Description	Mouse monoclonal antibody against Haemophilus Influenza Type b (Hib) conjugated to HRP
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Inactivated Haemophilus Influenza Type B
Preservatives	0.1M PBS, PH 8.0, with protein stabilizer
Storage	Store at 4°C for 6 months. It is recommended to aliquote to ~10ul per vial and store at -20°C for long term storage. Avoid repeated freezing and thawing cycles.
Note	For research use only
Expiration Date	12 months from date of receipt.