

Product Datasheet

Human CSF2RA protein (orb594713)

| | |
|--------------------------|---|
| Description | Recombinant Human Granulocyte-macrophage colony-stimulating factor receptor subunit alpha(CSF2RA),pa... |
| Endotoxins | Less than 1.0 EU/μg as determined by LAL method. |
| Preservatives | Lyophilized from a 0.2 μm filtered 20 mM PB, 150 mM NaCl, pH 7.4 |
| Form/Appearance | Lyophilized powder |
| Storage | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |
| Tag | C-terminal 6xHis-tagged |
| Note | For research use only |
| Application notes | Partial |
| Protein Sequence | EKSDLRTVAPASSLNVRFDSTRMNLWDCQENTTFSKCFLTDKKNRVVEPRLSN NECSCTFREICLHEGVTFEVHVNTSQRGFQQKLLYPNSGREGTAAQNFSCFIYNA DLMNCTWARGPTAPRDVQYFLYIRNSKRRREIRCPYYIQDSGTHVGCHLDNLG LTSRNYFLVNGTSREIGIQFFDSDLDTKKIERFNPPSNVTVRCNTHCLVRWKQPR TYQKLSYLDYQYQLDVHRKNTQPGTENLLINVSGDLENRYNFPSEPRAKHSVKI RAADVRLNWSSWSEAIEFGSDDG |
| Purity | Greater than 95% as determined by SDS-PAGE. |
| Source | Mammalian cell |
| MW | 35.5 kDa |
| Uniprot ID | P15509 |
| Expiration Date | 6 months from date of receipt. |