

Product Datasheet

Human FGFR3 protein (orb594741)

Description	Recombinant Human Fibroblast growth factor receptor 3(FGFR3),partial (Active)
Endotoxins	Less than 1.0 EU/μg as determined by LAL method.
Preservatives	Lyophilized from a 0.2 μm filtered 1xPBS, pH 7.4
Form/Appearance	Lyophilized powder
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Tag	C-terminal Fc-tagged
Note	For research use only
Application notes	Partial
Protein Sequence	ESLGTEQRVVGRAAEVPGPEPGQEQQLVFGSGDAVELSCPPPGGGPMGPTVW VKDGTGLVPSERVLVGPQRLQVLNASHEDSGAYSCRQRLTQRVLCHF SVRVTD APSSGDDDEGDEAEEDTGVDTGAPYWTRPERMDKLLAVPAANTVRFRC PAA GNPTPSISWLKNGREFRGEHRIGGIKLRHQWVSLVMESVVPDRGN YTCVVEN KFGSIRQTYTLDVLESPHRPILQAGLPANQTAVLGSDVEFHCKVYSDAQPHIQW LKHVEVNGSKVGPDPYVTVLKTAGANTTDKELEVELSLHNVT FEDAGEYTC LA GNSIGFSHHSAWLVVLPAAEELVEADEAGSVYAG
Purity	Greater than 95% as determined by SDS-PAGE.
Source	Mammalian cell
MW	64.8 kDa
Uniprot ID	P22607
Expiration Date	6 months from date of receipt.