

Product Datasheet

Human TGFB2 protein (orb594750)

Description	Recombinant Human Transforming growth factor beta-2(TGFB2) (Active)
Endotoxins	Less than 1.0 EU/μg as determined by LAL method.
Preservatives	Lyophilized from a 0.2 μm filtered 4 mM HCl
Form/Appearance	Lyophilized powder
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Tag	Tag-Free
Note	For research use only
Application notes	Full Length of Mature Protein
Protein Sequence	ALDAAYCFRNVQDNCCLRPLYIDFKRDLGWKWIHEPKGYNANFCAGACPYLWS SDTQHRSRVLSLYNTINPEASASPCCVSQDLEPLTILYIYGKTPKIEQLSNMIVKSK CS
Purity	Greater than 95% as determined by SDS-PAGE.
Source	Mammalian cell
MW	12.7 kDa
Uniprot ID	P61812
Expiration Date	6 months from date of receipt.