

Product Datasheet

Human FGF8 protein (orb594757)

Description	Recombinant Human Fibroblast growth factor 8(FGF8),partial (Active)
Endotoxins	Less than 1.0 EU/μg as determined by LAL method.
Preservatives	Lyophilized from a 0.2 μm filtered 1xPBS, pH 7.4
Form/Appearance	Lyophilized powder
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Tag	Tag-Free
Note	For research use only
Application notes	Partial of Isoform 3
Protein Sequence	QVTVQSSPNFTQHVREQSLVTDQLSRRLIRTYQLYSRTSGKHKVQVLANKRINAM AEDGDPFAKLIVETDTFGSRVRVARGAETGLYICMNKKGKLIKSNKGKDCVFTE IVLENNYTALQNAKYEGWYMAFTRKGRPRKGSKTRQHQREVHFMKRLPRGHHT TEQSLRFEFLNYPPFTRSLRGSQRTWAPEPR
Purity	Greater than 95% as determined by SDS-PAGE.
Source	E.coli
MW	22.5 kDa
Uniprot ID	P55075
Expiration Date	6 months from date of receipt.