

Product Datasheet

Mouse Kitlg protein (orb594761)

Description	Recombinant Mouse Kit ligand(Kitlg),partial (Active)
Endotoxins	Less than 1.0 EU/μg as determined by LAL method.
Preservatives	Lyophilized from a 0.2 μm filtered 1xPBS, pH 7.4
Form/Appearance	Lyophilized powder
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Tag	Tag-Free
Note	For research use only
Application notes	Partial
Protein Sequence	KEICGNPVTDNVKDITKLVANLPNDYMITLNYVAGMDVLP SHCWLRDMMVIQLSLS LTTLLDKFSNISEGLSNYSIIDKLGKIVDDLVL CMEENAPKNIKESPKRPETRSFTPE EFFSIFNRSIDAFKDFMVASDTSDCVLSSTLGPEKDSRVSVTKPFMLPPVA
Purity	Greater than 95% as determined by SDS-PAGE.
Source	E.coli
MW	18.4 kDa
Uniprot ID	P20826
Expiration Date	6 months from date of receipt.