

Product Datasheet

Mouse Ngf protein (orb594762)

Description	Recombinant Mouse Beta-nerve growth factor(Ngf),partial (Active)
Endotoxins	Less than 1.0 EU/μg as determined by LAL method.
Preservatives	Lyophilized from a 0.2 μm filtered solution of 20mM PB, 200mM NaCl, pH 8.0.
Form/Appearance	Lyophilized powder
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Tag	Tag-Free
Note	For research use only
Application notes	Partial
Protein Sequence	MGEFSVCDSVSVVWVGDKTTATDIKIGKEVTVLAEVNINNSVFRQYFFETKCRASN PVESGCRGIDSKHWNSYCTTTHTFVKALTTDEKQAAWRFIRIDTACVCLSRKAT R
Purity	Greater than 95% as determined by SDS-PAGE.
Source	E.coli
MW	12.4 kDa
Uniprot ID	P01139
Expiration Date	6 months from date of receipt.