

## Product Datasheet

### Mouse Fgf9 protein (orb594767)

<b>Description</b>	Recombinant Mouse Fibroblast growth factor 9(Fgf9) (Active)
<b>Endotoxins</b>	Less than 1.0 EU/μg as determined by LAL method.
<b>Preservatives</b>	0.2 μm filtered 20 mM Tris-HCl, 150 mM NaCl, 5% Trehalose, 1 mM EDTA, 20% Glycerol, 1 mM DTT, pH 8.5
<b>Form/Appearance</b>	Liquid
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Tag</b>	N-terminal 6xHis-tagged
<b>Note</b>	For research use only
<b>Application notes</b>	Full Length
<b>Protein Sequence</b>	MAPLGEVGSYFGVQDAVPFGNVPVLPVDSPLVLLSDHLGQSEAGGLPRGPAVTD LDHLKGILRRRQLYCRITGFHLEIFPNGTIQGRKDHSRFGILEFISIAVGLVSIRGV DSGLYLG MNEKGELYGSEKLTQECVFREQFEENWYNTYSSNLYKHVDTGRRYY VALNKDGTREGTRTKRHQKFTHFLPRPVPDPKVPPELYKDILSQS
<b>Purity</b>	Greater than 95% as determined by SDS-PAGE.
<b>Source</b>	E.coli
<b>MW</b>	24.4 kDa
<b>Uniprot ID</b>	<b>P54130</b>
<b>Expiration Date</b>	6 months from date of receipt.