

Product Datasheet

Human IL17RA protein (orb594789)

Description	Recombinant Human Interleukin-17 receptor A(IL17RA),partial (Active)
Endotoxins	Less than 1.0 EU/μg as determined by LAL method.
Preservatives	Lyophilized from a 0.2 μm filtered 1xPBS, pH 7.4
Form/Appearance	Lyophilized powder
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Tag	C-terminal Fc-tagged
Note	For research use only
Application notes	Extracellular Domain
Protein Sequence	LRLLDHRALVCSQPGLNCTVKNSTCLDDSWIHPRNLTSPSPKDLQQLHFAHTQ QGDLFPVAHIEWTLQTDASILYLEGAELSVLQLNTNERLCVRFEFLSKLRHHHRR WRFTFSHFVDPDQEYEVTVHHLPKPIPDGDPNHQSKNFLVPDCEHARMKVTTTP CMSSGSLWDPNITVETLEAHQLRVSFTLWNESTHYQILLTSFPHMENHSCFEHM HHIPAPRPEEFHQRSNVTLTRNLKGCCRHQVQIQPFSSCLNDCLRHSATVSCP EMPDTPEPIPDYMPPLW
Purity	Greater than 90% as determined by SDS-PAGE.
Source	Mammalian cell
MW	60.6 kDa
Uniprot ID	Q96F46
Expiration Date	6 months from date of receipt.