


## Product Datasheet

### Human IL2 protein (orb594792)

<b>Description</b>	Recombinant Human Interleukin-2(IL2) (Active)
<b>Endotoxins</b>	Less than 1.0 EU/μg as determined by LAL method.
<b>Preservatives</b>	Lyophilized from a 0.2 μm filtered 10 mM sodium citrate, pH 4.0
<b>Form/Appearance</b>	Lyophilized powder
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Tag</b>	Tag-Free
<b>Note</b>	For research use only
<b>Application notes</b>	Full Length of Mature Protein
<b>Protein Sequence</b>	APTSSSTKKTQLQLEHLLLDLQMLNGINNYKNPKLTRMLTFKFYMPKKATELKHL QCLEEELKPLEEVLNLAQSKNFHLRPRDLISINIVIVLELKGSETTFMCEYADETAT IVEFLNRWITFCQSIISTLT
<b>Purity</b>	Greater than 95% as determined by SDS-PAGE.
<b>Source</b>	E.coli
<b>MW</b>	15.5 kDa
<b>Uniprot ID</b>	<b>P60568</b>
<b>Expiration Date</b>	6 months from date of receipt.



(Tris-Glycine gel)  
Discontinuous  
SDS-PAG...