

Product Datasheet

Human FCGR3A protein (orb594865)

| | |
|--------------------------|---|
| Description | Recombinant Human Low affinity immunoglobulin gamma Fc region receptor III-A(FCGR3A),partial (Active... |
| Endotoxins | Less than 1.0 EU/μg as determined by LAL method. |
| Preservatives | Lyophilized from a 0.2 μm filtered 20 mM PB, 150 mM NaCl, pH 7.2 |
| Form/Appearance | Lyophilized powder |
| Storage | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |
| Tag | C-terminal 6xHis-tagged |
| Note | For research use only |
| Application notes | Extracellular Domain |
| Protein Sequence | GMRTEDLPKAVVFLEPQWYRVLEKDSVTLKQCQAYSPEDNSTQWFHNESSLISQ ASSYFIDAATVDDSGEYRCQTNLSTLSDPVQLEVHIGWLLLQAPRWVFKEEDPIH LRCHSWKNTALHKVTYLQNGKGRKYFHHNSDFYIPKATLKDSGSYFCRGLFGSK NVSSETVNITITQGLAVSTISSFFPPGYQ |
| Purity | Greater than 95% as determined by SDS-PAGE. |
| Source | Mammalian cell |
| MW | 22.61 kDa |
| Uniprot ID | P08637 |
| Expiration Date | 6 months from date of receipt. |