

Product Datasheet

Human SIRPA protein (orb594881)

Description	Recombinant Human Tyrosine-protein phosphatase non-receptor type substrate 1(SIRPA),partial (Active)
Endotoxins	Less than 1.0 EU/μg as determined by LAL method.
Preservatives	Lyophilized from a 0.2 μm filtered 20 mM PB, 150 mM NaCl, pH 7.2
Form/Appearance	Lyophilized powder
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Tag	C-terminal 6xHis-tagged
Note	For research use only
Application notes	Partial
Protein Sequence	EEELQVIQPKSVLVAAGETATLRCTATSLIPVGP IQWFRGAGPGRELIYNQKEG HFPRVTTVSDLTKRNNMDF SIRIGNITPADAGTYCVKFRKGPDDVEFKSGAGT ELSVRAKPSAPVVSGPAARATPQHTVSFTCESHGFSPRDITLKWFKNGNELSDF QTNVDPVGESVSYSIHSTAKVVLTR EDVHSQVICEVAHVTLQGDPLRGTANLSET IRVPPTLEVTQPVRAENQVNVTCQVRKFYPQRLQLTWLENGNVSRTETASTVT ENKDGTYNWM SWLLVNVSAHRDDVKLTCQVEHDGQPAVSKSHDLKVS AHPKE QGSNTAAENTGSNER
Purity	Greater than 95% as determined by SDS-PAGE.
Source	Mammalian cell
MW	38.31 kDa
Uniprot ID	P78324
Expiration Date	6 months from date of receipt.