

## Product Datasheet

### Virus mpt64 protein (orb595064)

<b>Description</b>	Recombinant Mycobacterium tuberculosis Immunogenic protein MPT64(mpt64)
<b>Reactivity</b>	Virus
<b>Target</b>	mpt64
<b>Preservatives</b>	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.
<b>Form/Appearance</b>	Liquid or Lyophilized powder
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Tag</b>	N-terminal 10xHis-tagged
<b>Note</b>	For research use only
<b>Application notes</b>	Full Length of Mature Protein
<b>Protein Sequence</b>	APKTYCEELKGTDTGQACQIQMSDPAYNINISLPSYYPDQKSLNIAQTRDKFLS AATSSTPREAPYELNITSATYQSAIPPRGTQAVVLKVYQNAGGTHPTTTYKAFDW DQAYRKPITYDTLWQADTDLPVVFPIVQGELSKQTGQVSIAPNAGLDPVNYQ NFAVTNDGVIFFFNPGELLPEAAGPTQVLVPRSAIDSMILA
<b>Purity</b>	Greater than 85% as determined by SDS-PAGE.
<b>Source</b>	Baculovirus
<b>MW</b>	24.9 kDa
<b>Uniprot ID</b>	<b>P9WIN9</b>
<b>Expiration Date</b>	6 months from date of receipt.