


Product Datasheet

PLC- γ 2 (phospho-Tyr1217) Antibody (orb6042)

Description	PLC- γ 2 (phospho-Tyr1217) Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	P16885
Expiration Date	12 months from date of receipt.



Immunohistochemical analysis of formalin...