

Product Datasheet

Human CA9 protein (orb604891)

Description	Recombinant Human Carbonic anhydrase 9(CA9),partial
Endotoxins	Not test.
Preservatives	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.
Form/Appearance	Liquid or Lyophilized powder
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Tag	N-terminal GST-tagged
Note	For research use only
Application notes	Partial
Protein Sequence	DLPSEEDSPREEDPPGEEDLPGEEDLPGEEDLPEVKPKSEEEGSLKLEDLPTVEA PGDPQEPQNNHRDKEGDDQSHWRYGGDPPWPRVSPACAGRFQSPVDIRPQL AAFCPALRPLELLGFQLPPLPELRLRNNNGHSVQLTLPPGLEMALGPGREYRALQL HLHWGAAGRPGSEHTVEGHRFP AEIHVVHLSTAFARVDEALGRPGLAVLAAFL EEGPEENSAYEQLLSRLEEIAEEGSETQVPGLDISALLPSDFSRYFYEGSLTTPPC AQGVIWTVFNQTVMLSAKQLHTLSDTLWGP GDSRLQLNFRATQPLNGRVIEASF PAGVDSSPR
Purity	Greater than 90% as determined by SDS-PAGE.
Source	E.coli
MW	63.8 kDa
Uniprot ID	Q16790
Expiration Date	6 months from date of receipt.