

Product Datasheet

Mouse Hadhb protein (orb605096)

Description	Recombinant Mouse Trifunctional enzyme subunit beta, mitochondrial(Hadhb)
Endotoxins	Not test.
Preservatives	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.
Form/Appearance	Liquid or Lyophilized powder
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Tag	N-terminal 10xHis-tagged and C-terminal Myc-tagged
Note	For research use only
Application notes	Full Length of Mature Protein
Protein Sequence	SAPAVQTKSKKTLAKPNMKNIVVVEGVRIPLLSGTSYKDLMPHDLARAALSGLL HRTNIPKDVVDYIIFGTVIQEVKTSNVAREEALGAGFSDKTPAHTVTMACISSNQA MTTAVGLIASGQCDVVVAGGVELMSDVPIRHSRNMKMLDLNKAKTLGQRLS LLSKFRLNFLSPELPAVAEFSTNETMGHSADRLAAAFVSRMEQDEYALRSHSLA KKAQDEGHLSDIVPFKVPKGDVTVDKNGIRPSSLEQMAKLPKPAFIKPYGTVTAAN SSFLTGDGASAMLIMSEDRALAMGYKPKAYLRDFIYVSQDPKDQLLLGPTYATPKV LEKAGLTMNDIDAFEFHEAFSGQILANFKAMDSDFWFAQNYMGRKTKVGSPPLEK FNIWGGSLSLGHFPFGATGCRMLVMAAANRLRKDGGQYALVAACAAGGQGHAMI VEAYPK
Purity	Greater than 85% as determined by SDS-PAGE.
Source	E.coli
MW	54.6 kDa
Uniprot ID	Q99JY0
Expiration Date	6 months from date of receipt.