

Product Datasheet

Mouse Properdin protein (orb605204)

Description	Recombinant Mouse Properdin(Cfp)
Endotoxins	Not test.
Preservatives	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.
Form/Appearance	Liquid or Lyophilized powder
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Tag	N-terminal 10xHis-tagged and C-terminal Myc-tagged
Note	For research use only
Application notes	Full Length of Mature Protein
Protein Sequence	SDPVLCTFYEESSGRCKGLLGRDIRVEDCCLNAAAFQEHGGLCQACRSPQ WSAWSLWGPCSVTCSEGSQLRHRRCVGRGGQCSENVAPGTLEWQLQACEDQ PCCPEMGGWSEWGPWGPCSVTCSKGTQIRQRVCDNPAPKCGGHCPGEAQQS QACDTQKTCPTHGAWASWGPWSPCSGSLGGAQEPKETRSRSCSAPAPSHQP PGKPCSGPAYEHKACSGLPVAGGWGPWSPLSPCSVTCGLGQTLEQRTCDH PAPRHGGPFCAGDATRNQMCNKAVPCPVNGEWEAWGKWSDCSRLRMSINCE GTPGQQSRSRSCGGRKFNGKPCAGKLQDIRHCYNHNCIMKGSWSQWSTWSL CTPPCSPNATRVQRQLCTPLLPKYPTVSMVEGQGEKNVTFWGTTPRPLCEALQG QKLVEEKRSLHVPVCKDPEEKKP
Purity	Greater than 85% as determined by SDS-PAGE.
Source	Mammalian cell
MW	53.0 kDa
Uniprot ID	P11680
Expiration Date	6 months from date of receipt.