

Product Datasheet

Human MADCAM1 protein (orb605268)

Description	Recombinant Human Mucosal addressin cell adhesion molecule 1(MADCAM1),partial
Endotoxins	Not test.
Preservatives	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.
Form/Appearance	Liquid or Lyophilized powder
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Tag	N-terminal 10xHis-tagged
Note	For research use only
Application notes	Extracellular Domain
Protein Sequence	QSLQVKPLQVEPPEPVVAVALGASRQLTCRLACADRGASVQWRGLDTSLGAVQ SDTGRSVLTVRNASLSAAGTRVCVCGSCGGRTFQHTVQLLVYAFPDQLTVSPAAL VPGDPEVACTAHKVTPVDPNALSFSLLVGGQELEGAQALGPEVQEEEEEPQGDE DVLFRVTERWRLPPLGTPVPALYQCQATMRLPGLLESHRQAIPVLHSPTSPEPPD TTSPEPDTTSPESPDTTSPQPPDTSPEPPDKTSPEPAPQQGSTHTPRSPGSTR TRRPEISQAGPTQGEVIPTGSSKPADGQ
Purity	Greater than 85% as determined by SDS-PAGE.
Source	Mammalian cell
MW	34.9 kDa
Uniprot ID	Q13477
Expiration Date	6 months from date of receipt.